

HOTBORN

UT-500X

UROLOGICAL TABLE(Lithotripsy table)

Table compatible with Extracorporeal Shockwave Lithotripter for Lithotripsy surgery.



FEATURE

- ▶ Micro-motors and accurate slide ensure quick, smooth, and balance table movement.
- ▶ Carbon fibre table top allows clinicians to get high-quality, clear imaging with low X-ray doses.
- ▶ Five movement table: height, lateral tilt, trendelenburg, longitudinal, lateral travel.
- ▶ The ESWL section is special for Extracorporeal Shockwave Lithotripter.
- ▶ The base offers improved ergonomics and surgeon access.
- ▶ High-quality big casters and braking system allow clinicians to easily move the table.
- ▶ 304 stainless steel design for durability and easy cleaning
- ▶ Optional battery which can work one week (60~80 operations)
- ▶ Optional function level position.
- ▶ High-quality material, advanced and mature technology ensure the quality.



SPECIFICATIONS

- ▶ Table top(L*W): 2500*620 mm
- ▶ ESWL section: 480*280 (arc shaped)
- ▶ Longitudinal travel: 210mm
- ▶ Lateral travel: 160mm
- ▶ Height: 800~1080 mm
- ▶ Trendelenburg& Reverse~: 19°
- ▶ Lateral tilt: 16°
- ▶ Load capacity: 250 kg
- ▶ Power: 110V~220V 50Hz, 350VA



STANDARD ACCESSORIES

- ▶ P-001 Mattress- 1 set
- ▶ P-002 Cable controller-1 set
- ▶ P-003 Shoulder & Lateral support-1 pair
- ▶ P-005 Arm board- 1 pair
- ▶ P-006 Arm fixed belt-1 pair
- ▶ P-007 Fixed belt-1 pair
- ▶ P-008 IV Pole- 1 set
- ▶ P-009 Anaesthesia screen-1 set

OPTIONAL ACCESSORIES

- ▶ P-010 Leg support
- ▶ P-014 Hanging arm board
- ▶ P-015 Double arm board